

## Antibodies Product Datasheet

|                               |                                                                                                                                                                      |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | CD84 Antibody                                                                                                                                                        |
| <b>Catalog Number:</b>        | C30386                                                                                                                                                               |
| <b>Reactivity:</b>            | Human, Mouse                                                                                                                                                         |
| <b>Applications:</b>          | WB, IHC-p, IF(paraffin section), ELISA                                                                                                                               |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.                                                                 |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                                     |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                    |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from the Internal region of human CD84. AA range:131-180                                              |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                              |
| <b>Gene Name:</b>             | CD84                                                                                                                                                                 |
| <b>Protein Name:</b>          | SLAM family member 5                                                                                                                                                 |
| <b>Other Name:</b>            | CD84; SLAMF5; SLAM family member 5; Cell surface antigen MAX.3; Hly9-beta; Leukocyte differentiation antigen CD84; Signaling lymphocytic activation molecule 5; CD84 |
| <b>Human Swiss Prot No.:</b>  | Q9UIB8                                                                                                                                                               |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                              |
| <b>Specificity:</b>           | CD84 Polyclonal Antibody detects endogenous levels of CD84 protein.                                                                                                  |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                         |