

Antibodies Product Datasheet

Product Name:	CD85c Antibody
Catalog Number:	C30388
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human LILRB5. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	LILRB5
Protein Name:	Leukocyte immunoglobulin-like receptor subfamily B member 5
Other Name:	LILRB5; LIR8; Leukocyte immunoglobulin-like receptor subfamily B member 5; CD85 antigen-like family member C; Leukocyte immunoglobulin-like receptor 8; LIR-8; CD85c
Human Swiss Prot No.:	O75023
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD85c Polyclonal Antibody detects endogenous levels of CD85c protein.
Storage and Stability:	-20°C/1 year