

Antibodies Product Datasheet

Product Name:	CD85d Antibody
Catalog Number:	C30389
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human LILRB2. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	LILRB2
Protein Name:	Leukocyte immunoglobulin-like receptor subfamily B member 2
Other Name:	LILRB2; ILT4; LIR2; MIR10; Leukocyte immunoglobulin-like receptor subfamily B member 2; LIR-2; Leukocyte immunoglobulin-like receptor 2; CD85 antigen-like family member D; Immunoglobulin-like transcript 4; ILT-4; Monocyte/macrophage immunoglobulin-like re
Human Swiss Prot No.:	Q8N423
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD85d Polyclonal Antibody detects endogenous levels of CD85d protein.
Storage and Stability:	-20°C/1 year