

Antibodies Product Datasheet

Product Name:	CD85f Antibody
Catalog Number:	C30391
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human LILRA5. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	LILRA5
Protein Name:	Leukocyte immunoglobulin-like receptor subfamily A member 5
Other Name:	LILRA5; ILT11; LILRB7; LIR9; Leukocyte immunoglobulin-like receptor subfamily A member 5; CD85 antigen-like family member F; Immunoglobulin-like transcript 11; ILT-11; Leukocyte immunoglobulin-like receptor 9; LIR-9; CD85f
Human Swiss Prot No.:	A6NI73
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD85f Polyclonal Antibody detects endogenous levels of CD85f protein.
Storage and Stability:	-20°C/1 year