

Antibodies Product Datasheet

Product Name:	CD85g Antibody
Catalog Number:	C30392
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human LILRA4. AA range:391-440
Concentration:	1 mg/ml
Gene Name:	LILRA4
Protein Name:	Leukocyte immunoglobulin-like receptor subfamily A member 4
Other Name:	LILRA4; ILT7; Leukocyte immunoglobulin-like receptor subfamily A member 4; CD85 antigen-like family member G; Immunoglobulin-like transcript 7; ILT-7; CD85g
Human Swiss Prot No.:	P59901
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD85g Polyclonal Antibody detects endogenous levels of CD85g protein.
Storage and Stability:	-20°C/1 year