

Antibodies Product Datasheet

Product Name:	LIR-7 Antibody
Catalog Number:	C30393
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human LILRA2. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	LILRA2
Protein Name:	Leukocyte immunoglobulin-like receptor subfamily A member 2
Other Name:	LILRA2; ILT1; LIR7; Leukocyte immunoglobulin-like receptor subfamily A member 2; CD85 antigen-like family member H; Immunoglobulin-like transcript 1; ILT-1; Leukocyte immunoglobulin-like receptor 7; LIR-7; CD85h
Human Swiss Prot No.:	Q8N149
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LIR-7 Polyclonal Antibody detects endogenous levels of LIR-7 protein.
Storage and Stability:	-20°C/1 year