

## Antibodies Product Datasheet

<b>Product Name:</b>	LIR-1 Antibody
<b>Catalog Number:</b>	C30395
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from the N-terminal region of human LILRB1. AA range:21-70
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	LILRB1
<b>Protein Name:</b>	Leukocyte immunoglobulin-like receptor subfamily B member 1
<b>Other Name:</b>	LILRB1; ILT2; LIR1; MIR7; Leukocyte immunoglobulin-like receptor subfamily B member 1; LIR-1; Leukocyte immunoglobulin-like receptor 1; CD85 antigen-like family member J; Immunoglobulin-like transcript 2; ILT-2; Monocyte/macrophage immunoglobulin-like rec
<b>Human Swiss Prot No.:</b>	Q8NHL6
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	LIR-1 Polyclonal Antibody detects endogenous levels of LIR-1 protein.
<b>Storage and Stability:</b>	-20°C/1 year