

Antibodies Product Datasheet

Product Name:	CD93 Antibody
Catalog Number:	C30398
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD93. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	CD93
Protein Name:	Complement component C1q receptor
Other Name:	CD93; C1QR1; MXRA4; Complement component C1q receptor; C1q/MBL/SPA receptor; C1qR; C1qR(p); C1qRp; CDw93; Complement component 1 q subcomponent receptor 1; Matrix-remodeling-associated protein 4; CD93
Human Swiss Prot No.:	Q9NPY3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD93 Polyclonal Antibody detects endogenous levels of CD93 protein.
Storage and Stability:	-20°C/1 year