

Antibodies Product Datasheet

Product Name:	CD96 Antibody
Catalog Number:	C30401
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD96. AA range:291-340
Concentration:	1 mg/ml
Gene Name:	CD96
Protein Name:	T-cell surface protein tactile
Other Name:	CD96; T-cell surface protein tactile; Cell surface antigen CD96; T cell-activated increased late expression protein; CD96
Human Swiss Prot No.:	P40200
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD96 Polyclonal Antibody detects endogenous levels of CD96 protein.
Storage and Stability:	-20°C/1 year