

Antibodies Product Datasheet

Product Name:	CD98 Antibody
Catalog Number:	C30403
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human SLC3A2. AA range:491-540
Concentration:	1 mg/ml
Gene Name:	SLC3A2
Protein Name:	4F2 cell-surface antigen heavy chain
Other Name:	SLC3A2; MDU1; 4F2 cell-surface antigen heavy chain; 4F2hc; 4F2 heavy chain antigen; Lymphocyte activation antigen 4F2 large subunit; CD98
Human Swiss Prot No.:	P08195
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD98 Polyclonal Antibody detects endogenous levels of CD98 protein.
Storage and Stability:	-20°C/1 year