

Antibodies Product Datasheet

Product Name:	CD101 Antibody
Catalog Number:	C30406
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD101. AA range:731-780
Concentration:	1 mg/ml
Gene Name:	CD101
Protein Name:	Immunoglobulin superfamily member 2
Other Name:	CD101; EWI101; IGSF2; V7; Immunoglobulin superfamily member 2; IgSF2; Cell surface glycoprotein V7; Glu-Trp-Ile EWI motif-containing protein 101; EWI-101; CD101
Human Swiss Prot No.:	Q93033
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD101 Polyclonal Antibody detects endogenous levels of CD101 protein.
Storage and Stability:	-20°C/1 year