

Antibodies Product Datasheet

Product Name:	IL-1RII Antibody
Catalog Number:	C30423
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human IL1R2. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	IL1R2
Protein Name:	Interleukin-1 receptor type 2
Other Name:	IL1R2; IL1RB; Interleukin-1 receptor type 2; IL-1R-2; IL-1RT-2; IL-1RT2; CD121 antigen-like family member B; CDw121b; IL-1 type II receptor; Interleukin-1 receptor beta; IL-1R-beta; Interleukin-1 receptor type II; CD121b
Human Swiss Prot No.:	P27930
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-1RII Polyclonal Antibody detects endogenous levels of IL-1RII protein.
Storage and Stability:	-20°C/1 year