

Antibodies Product Datasheet

Product Name:	EMMPRIN Antibody
Catalog Number:	C30438
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human BSG. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	BSG
Protein Name:	Basigin
Other Name:	BSG; Basigin; 5F7; Collagenase stimulatory factor; Extracellular matrix metalloproteinase inducer; EMMPRIN; Leukocyte activation antigen M6; OK blood group antigen; Tumor cell-derived collagenase stimulatory factor; TCSF; CD147
Human Swiss Prot No.:	P35613
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EMMPRIN Polyclonal Antibody detects endogenous levels of EMMPRIN protein.
Storage and Stability:	-20°C/1 year