

Antibodies Product Datasheet

Product Name:	Blr1 Antibody
Catalog Number:	C30442
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human CXCR5. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	CXCR5
Protein Name:	C-X-C chemokine receptor type 5
Other Name:	CXCR5; BLR1; MDR15; C-X-C chemokine receptor type 5; CXC-R5; CXCR-5; Burkitt lymphoma receptor 1; Monocyte-derived receptor 15; MDR-15; CD185
Human Swiss Prot No.:	P32302
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Blr1 Polyclonal Antibody detects endogenous levels of Blr1 protein.
Storage and Stability:	-20°C/1 year