

Antibodies Product Datasheet

Product Name:	IL-18R Antibody
Catalog Number:	C30446
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human IL18R1. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	IL18R1
Protein Name:	Interleukin-18 receptor 1
Other Name:	IL18R1; IL1RRP; Interleukin-18 receptor 1; IL-18R-1; IL-18R1; CD218 antigen-like family member A; CDw218a; IL1 receptor-related protein; IL-1Rrp; IL1R-rp; CD218a
Human Swiss Prot No.:	Q13478
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-18R Polyclonal Antibody detects endogenous levels of IL-18R protein.
Storage and Stability:	-20°C/1 year