

Antibodies Product Datasheet

Product Name:	IL-12R β 1 Antibody
Catalog Number:	C30455
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human IL12RB1. AA range:211-260
Concentration:	1 mg/ml
Gene Name:	IL12RB1
Protein Name:	Interleukin-12 receptor subunit beta-1
Other Name:	IL12RB1; IL12R; IL12RB; Interleukin-12 receptor subunit beta-1; IL-12 receptor subunit beta-1; IL-12R subunit beta-1; IL-12R-beta-1; IL-12RB1; IL-12 receptor beta component; CD212
Human Swiss Prot No.:	P42701
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-12R β 1 Polyclonal Antibody detects endogenous levels of IL-12R β 1 protein.
Storage and Stability:	-20°C/1 year