

Antibodies Product Datasheet

Product Name:	CD244 Antibody
Catalog Number:	C30457
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD244. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	CD244
Protein Name:	Natural killer cell receptor 2B4
Other Name:	CD244; 2B4; Natural killer cell receptor 2B4; NK cell activation-inducing ligand; NAIL; NK cell type I receptor protein 2B4; NKR2B4; h2B4; CD244 Q9BZW8
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD244 Polyclonal Antibody detects endogenous levels of CD244 protein.
Storage and Stability:	-20°C/1 year