

Antibodies Product Datasheet

Product Name:	CD307 Antibody
Catalog Number:	C30458
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human FCRL5. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	FCRL5
Protein Name:	Fc receptor-like protein 5
Other Name:	FCRL5; FCRH5; IRTA2; Fc receptor-like protein 5; FcR-like protein 5; FcRL5; BXMAS1; Fc receptor homolog 5; FcRH5; Immune receptor translocation-associated protein 2; CD307e
Human Swiss Prot No.:	Q96RD9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD307 Polyclonal Antibody detects endogenous levels of CD307 protein.
Storage and Stability:	-20°C/1 year