

Antibodies Product Datasheet

Product Name:	IL-13R α 2 Antibody
Catalog Number:	C30462
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Interleukin-13 receptor subunit alpha-2 at AA range: 241-290
Concentration:	1 mg/ml
Gene Name:	IL13RA2
Protein Name:	Interleukin-13 receptor subunit alpha-2
Other Name:	IL13RA2; IL13R; Interleukin-13 receptor subunit alpha-2; IL-13 receptor subunit alpha-2; IL-13R subunit alpha-2; IL-13R-alpha-2; IL-13RA2; Interleukin-13-binding protein; CD213a2
Human Swiss Prot No.:	Q14627
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-13R α 2 Polyclonal Antibody detects endogenous levels of IL-13R α 2 protein.
Storage and Stability:	-20°C/1 year