

## Antibodies Product Datasheet

|                               |                                                                                                                                                                   |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Bonzo Antibody                                                                                                                                                    |
| <b>Catalog Number:</b>        | C30466                                                                                                                                                            |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                                                                 |
| <b>Applications:</b>          | WB, IHC-p, IF(paraffin section), ELISA                                                                                                                            |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.                                                           |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                                  |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                 |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                             |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from the N-terminal region of human CXCR6. AA range:1-50                                           |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                           |
| <b>Gene Name:</b>             | CXCR6                                                                                                                                                             |
| <b>Protein Name:</b>          | C-X-C chemokine receptor type 6                                                                                                                                   |
| <b>Other Name:</b>            | CXCR6; BONZO; STRL33; TYMSTR; C-X-C chemokine receptor type 6; CXC-R6; CXCR-6; CDw186; G-protein coupled receptor STRL33; G-protein coupled receptor bonzo; CD186 |
| <b>Human Swiss Prot No.:</b>  | O00574                                                                                                                                                            |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                           |
| <b>Specificity:</b>           | Bonzo Polyclonal Antibody detects endogenous levels of Bonzo protein.                                                                                             |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                      |