

Antibodies Product Datasheet

Product Name:	CD169 Antibody
Catalog Number:	C30468
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human SIGLEC1. AA range:1321-1370
Concentration:	1 mg/ml
Gene Name:	SIGLEC1
Protein Name:	Sialoadhesin
Other Name:	SIGLEC1; SN; Sialoadhesin; Sialic acid-binding Ig-like lectin 1; Siglec-1; CD169
Human Swiss Prot No.:	Q9BZZ2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD169 Polyclonal Antibody detects endogenous levels of CD169 protein.
Storage and Stability:	-20°C/1 year