

Antibodies Product Datasheet

Product Name:	CD158a Antibody
Catalog Number:	C30474
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human KIR2DL1. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	KIR2DL1
Protein Name:	Killer cell immunoglobulin-like receptor 2DL1
Other Name:	KIR2DL1; CD158A; NKAT1; Killer cell immunoglobulin-like receptor 2DL1; CD158 antigen-like family member A; MHC class I NK cell receptor; Natural killer-associated transcript 1; NKAT-1; p58 natural killer cell receptor clones CL-42/47.11; p58 NK receptor C
Human Swiss Prot No.:	P43626
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD158a Polyclonal Antibody detects endogenous levels of CD158a protein.
Storage and Stability:	-20°C/1 year