

Antibodies Product Datasheet

Product Name:	CD158k Antibody
Catalog Number:	C30475
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p 1:50-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 221-270
Concentration:	1 mg/ml
Gene Name:	KIR3DL2 CD158K NKAT4
Protein Name:	Killer cell immunoglobulin-like receptor 3DL2 (CD158 antigen-like family member K) (MHC class I NK cell receptor) (Natural killer-associated transcript 4) (NKAT-4) (p70 natural killer cell receptor clone CL-5) (p70 NK receptor CL-5) (CD antigen CD158k)
Other Name:	Killer cell immunoglobulin-like receptor 3DL2 (CD158 antigen-like family member K) (MHC class I NK cell receptor) (Natural killer-associated transcript 4) (NKAT-4) (p70 natural killer cell receptor clone CL-5) (p70 NK receptor CL-5) (CD antigen CD158k)
Human Swiss Prot No.:	P43630
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous CD158k
Storage and Stability:	-20°C/1 year