

Antibodies Product Datasheet

Product Name:	CD158f1/2 Antibody
Catalog Number:	C30476
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p 1:50-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 121-170
Concentration:	1 mg/ml
Gene Name:	KIR2DL5A CD158F CD158F1 KIR2DL5 KIR2DL5B CD158F CD158F2 KIR2DL5 KIR2DLX
Protein Name:	Killer cell immunoglobulin-like receptor 2DL5A/B (CD antigen CD158f1/2)
Other Name:	
Human Swiss Prot No.:	Q8N109/Q8NHK3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous CD158f1/2
Storage and Stability:	-20°C/1 year