

Antibodies Product Datasheet

Product Name:	CD300c Antibody
Catalog Number:	C30477
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD300C. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	CD300C
Protein Name:	CMRF35-like molecule 6
Other Name:	CD300C; CMRF35; CMRF35A; CMRF35A1; IGSF16; CMRF35-like molecule 6; CLM-6; CD300 antigen-like family member C; CMRF35-A1; CMRF-35; Immunoglobulin superfamily member 16; IgSF16; CD300c
Human Swiss Prot No.:	Q08708
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD300c Polyclonal Antibody detects endogenous levels of CD300c protein.
Storage and Stability:	-20°C/1 year