

Antibodies Product Datasheet

Product Name:	CD158z Antibody
Catalog Number:	C30485
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human KIR3DL3. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	KIR3DL3
Protein Name:	Killer cell immunoglobulin-like receptor 3DL3
Other Name:	KIR3DL3; CD158Z; KIR3DL7; KIRC1; Killer cell immunoglobulin-like receptor 3DL3; CD158 antigen-like family member Z; Killer cell inhibitory receptor 1; CD158z
Human Swiss Prot No.:	Q8N743
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD158z Polyclonal Antibody detects endogenous levels of CD158z protein.
Storage and Stability:	-20°C/1 year