

Antibodies Product Datasheet

Product Name:	DDR2 Antibody
Catalog Number:	C30492
Reactivity:	Human, Mouse
Applications:	IF/ICC, IHC-p, ELISA
Dilutions:	IHC-p 1:50-200, IF1: 500 ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 31-80
Concentration:	1 mg/ml
Gene Name:	DDR2 NTRKR3 TKT TYRO10
Protein Name:	Discoidin domain-containing receptor 2 (Discoidin domain receptor 2) (EC 2.7.10.1) (CD167 antigen-like family member B) (Discoidin domain-containing receptor tyrosine kinase 2) (Neurotrophic tyrosine kinase, receptor-related 3) (Receptor protein-tyrosine
Other Name:	Discoidin domain-containing receptor 2 (Discoidin domain receptor 2) (EC 2.7.10.1) (CD167 antigen-like family member B) (Discoidin domain-containing receptor tyrosine kinase 2) (Neurotrophic tyrosine kinase, receptor-related 3) (Receptor protein-tyrosine
Human Swiss Prot No.:	Q16832
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous DDR2
Storage and Stability:	-20°C/1 year