

Antibodies Product Datasheet

Product Name:	CD218b Antibody
Catalog Number:	C30495
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p 1:50-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 291-340
Concentration:	1 mg/ml
Gene Name:	IL18RAP IL1R7
Protein Name:	Interleukin-18 receptor accessory protein (IL-18 receptor accessory protein) (IL-18RAcP) (Accessory protein-like) (AcPL) (CD218 antigen-like family member B) (CDw218b) (IL-1R accessory protein-like) (IL-1RAcPL) (Interleukin-1 receptor 7) (IL-1R-7) (IL-1R7)
Other Name:	Interleukin-18 receptor accessory protein (IL-18 receptor accessory protein) (IL-18RAcP) (Accessory protein-like) (AcPL) (CD218 antigen-like family member B) (CDw218b) (IL-1R accessory protein-like) (IL-1RAcPL) (Interleukin-1 receptor 7) (IL-1R-7) (IL-1R7)
Human Swiss Prot No.:	O95256
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous CD218b
Storage and Stability:	-20°C/1 year