

Antibodies Product Datasheet

Product Name:	CD298 Antibody
Catalog Number:	C30501
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human ATP1B3. AA range:222-271
Concentration:	1 mg/ml
Gene Name:	ATP1B3
Protein Name:	Sodium/potassium-transporting ATPase subunit beta-3
Other Name:	ATP1B3; Sodium/potassium-transporting ATPase subunit beta-3; Sodium/potassium-dependent ATPase subunit beta-3; ATPB-3; CD298
Human Swiss Prot No.:	P54709
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD298 Polyclonal Antibody detects endogenous levels of CD298 protein.
Storage and Stability:	-20°C/1 year