

Antibodies Product Datasheet

Product Name:	CD329 Antibody
Catalog Number:	C30515
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human SIGLEC8. AA range:81-130
Concentration:	1 mg/ml
Gene Name:	SIGLEC8
Protein Name:	Sialic acid-binding Ig-like lectin 8
Other Name:	SIGLEC8; SAF2; Sialic acid-binding Ig-like lectin 8; Siglec-8; CDw329; Sialoadhesin family member 2; SAF-2; CD329
Human Swiss Prot No.:	Q9NYZ4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD329 Polyclonal Antibody detects endogenous levels of CD329 protein.
Storage and Stability:	-20°C/1 year