

Antibodies Product Datasheet

Product Name:	CD239 Antibody
Catalog Number:	C30518
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human BCAM. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	BCAM
Protein Name:	Basal cell adhesion molecule
Other Name:	BCAM; LU; MSK19; Basal cell adhesion molecule; Auberger B antigen; B-CAM cell surface glycoprotein; F8/G253 antigen; Lutheran antigen; Lutheran blood group glycoprotein; CD239
Human Swiss Prot No.:	P50895
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD239 Polyclonal Antibody detects endogenous levels of CD239 protein.
Storage and Stability:	-20°C/1 year