

Antibodies Product Datasheet

Product Name:	CD225 Antibody
Catalog Number:	C30521
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human IFITM1. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	IFITM1
Protein Name:	Interferon-induced transmembrane protein 1
Other Name:	IFITM1; CD225; IFI17; Interferon-induced transmembrane protein 1; Dispanin subfamily A member 2a; DSPA2a; Interferon-induced protein 17; Interferon-inducible protein 9-27; Leu-13 antigen; CD225
Human Swiss Prot No.:	P13164
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD225 Polyclonal Antibody detects endogenous levels of CD225 protein.
Storage and Stability:	-20°C/1 year