

## Antibodies Product Datasheet

|                               |                                                                                                                                                          |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | IL-8RA Antibody                                                                                                                                          |
| <b>Catalog Number:</b>        | C30526                                                                                                                                                   |
| <b>Reactivity:</b>            | Human                                                                                                                                                    |
| <b>Applications:</b>          | WB, IHC-p, IF(paraffin section), ELISA                                                                                                                   |
| <b>Dilutions:</b>             | IHC-p 1:50-200, ELISA 1:10000-20000                                                                                                                      |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                         |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                        |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                    |
| <b>Immunogen:</b>             | Synthetic peptide from human protein at AA range: 1-50                                                                                                   |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                  |
| <b>Gene Name:</b>             | CXCR1 CMKAR1 IL8RA                                                                                                                                       |
| <b>Protein Name:</b>          | C-X-C chemokine receptor type 1 (CXC-R1) (CXCR-1) (CDw128a) (High affinity interleukin-8 receptor A) (IL-8R A) (IL-8 receptor type 1) (CD antigen CD181) |
| <b>Other Name:</b>            | C-X-C chemokine receptor type 1 (CXC-R1) (CXCR-1) (CDw128a) (High affinity interleukin-8 receptor A) (IL-8R A) (IL-8 receptor type 1) (CD antigen CD181) |
| <b>Human Swiss Prot No.:</b>  | P25024                                                                                                                                                   |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                  |
| <b>Specificity:</b>           | The antibody detects endogenous IL-8RA                                                                                                                   |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                             |