

Antibodies Product Datasheet

Product Name:	CD328 Antibody
Catalog Number:	C30534
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human SIGLEC7. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	SIGLEC7
Protein Name:	Sialic acid-binding Ig-like lectin 7
Other Name:	SIGLEC7; AIRM1; Sialic acid-binding Ig-like lectin 7; Siglec-7; Adhesion inhibitory receptor molecule 1; AIRM-1; CDw328; D-siglec; QA79 membrane protein; p75; CD328
Human Swiss Prot No.:	Q9Y286
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD328 Polyclonal Antibody detects endogenous levels of CD328 protein.
Storage and Stability:	-20°C/1 year