

Antibodies Product Datasheet

Product Name:	CD238 Antibody
Catalog Number:	C30538
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human KEL. AA range:351-400
Concentration:	1 mg/ml
Gene Name:	KEL
Protein Name:	CD238
Other Name:	Kell blood group, metallo-endopeptidase
Human Swiss Prot No.:	P23276
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD238 Polyclonal Antibody detects endogenous levels of CD238
Storage and Stability:	-20°C/1 year