

Antibodies Product Datasheet

Product Name:	CD209 Antibody
Catalog Number:	C30540
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from CD209 antigen at AA range: 261-310
Concentration:	1 mg/ml
Gene Name:	CD209
Protein Name:	CD209 antigen
Other Name:	CD209; CLEC4L; CD209 antigen; C-type lectin domain family 4 member L; Dendritic cell-specific ICAM-3-grabbing non-integrin 1; DC-SIGN; DC-SIGN1; CD209
Human Swiss Prot No.:	Q9NNX6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD209 Polyclonal Antibody detects endogenous levels of CD209 protein.
Storage and Stability:	-20°C/1 year