

Antibodies Product Datasheet

Product Name:	CD305 Antibody
Catalog Number:	C30541
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human LAIR1. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	LAIR1
Protein Name:	Leukocyte-associated immunoglobulin-like receptor 1
Other Name:	LAIR1; CD305; Leukocyte-associated immunoglobulin-like receptor 1; LAIR-1; hLAIR1; CD305
Human Swiss Prot No.:	Q6GTX8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD305 Polyclonal Antibody detects endogenous levels of CD305 protein.
Storage and Stability:	-20°C/1 year