

Antibodies Product Datasheet

Product Name:	CD177 Antibody
Catalog Number:	C30544
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD177. AA range:361-410
Concentration:	1 mg/ml
Gene Name:	CD177
Protein Name:	CD177 antigen
Other Name:	CD177; NB1; PRV1; CD177 antigen; Human neutrophil alloantigen 2a; HNA-2a; NB1 glycoprotein; NB1 GP; Polycythemia rubra vera protein 1; PRV-1; CD177
Human Swiss Prot No.:	Q8N6Q3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD177 Polyclonal Antibody detects endogenous levels of CD177 protein.
Storage and Stability:	-20°C/1 year