

Antibodies Product Datasheet

Product Name:	CD163b Antibody
Catalog Number:	C30545
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD163L1. AA range:391-440
Concentration:	1 mg/ml
Gene Name:	CD163L1
Protein Name:	Scavenger receptor cysteine-rich type 1 protein M160
Other Name:	CD163L1; CD163B; M160; Scavenger receptor cysteine-rich type 1 protein M160; CD163 antigen-like 1; CD163b
Human Swiss Prot No.:	Q9NR16
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD163b Polyclonal Antibody detects endogenous levels of CD163b protein.
Storage and Stability:	-20°C/1 year