

Antibodies Product Datasheet

Product Name:	CD276 Antibody
Catalog Number:	C30551
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD276. AA range:271-320
Concentration:	1 mg/ml
Gene Name:	CD276
Protein Name:	CD276 antigen
Other Name:	CD276; B7H3; PSEC0249; CD276 antigen; 4Ig-B7-H3; B7 homolog 3; B7-H3; Costimulatory molecule; CD276
Human Swiss Prot No.:	Q5ZPR3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD276 Polyclonal Antibody detects endogenous levels of CD276 protein.
Storage and Stability:	-20°C/1 year