

Antibodies Product Datasheet

Product Name:	BTN3A1/2/3 Antibody
Catalog Number:	C30554
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p 1:50-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 41-90
Concentration:	1 mg/ml
Gene Name:	BTN3A1/2/3 BTF3/4/5
Protein Name:	Butyrophilin subfamily 3 member A1/2/3 (CD antigen CD277)
Other Name:	
Human Swiss Prot No.:	O00481/P78410/O00478
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous BTN3A1/2/3
Storage and Stability:	-20°C/1 year