

Antibodies Product Datasheet

Product Name:	CD315 Antibody
Catalog Number:	C30561
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human PTGFRN. AA range:601-650
Concentration:	1 mg/ml
Gene Name:	PTGFRN
Protein Name:	Prostaglandin F2 receptor negative regulator
Other Name:	PTGFRN; CD9P1; EWIF; FPRP; KIAA1436; Prostaglandin F2 receptor negative regulator; CD9 partner 1; CD9P-1; Glu-Trp-Ile EWI motif-containing protein F; EWI-F; Prostaglandin F2-alpha receptor regulatory protein; Prostaglandin F2-alpha receptor-associated pro
Human Swiss Prot No.:	Q9P2B2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD315 Polyclonal Antibody detects endogenous levels of CD315 protein.
Storage and Stability:	-20°C/1 year