

Antibodies Product Datasheet

Product Name:	CD316 Antibody
Catalog Number:	C30567
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human IGSF8. AA range:451-500
Concentration:	1 mg/ml
Gene Name:	IGSF8
Protein Name:	Immunoglobulin superfamily member 8
Other Name:	IGSF8; CD81P3; EWI2; KCT4; Immunoglobulin superfamily member 8; IgSF8; CD81 partner 3; Glu-Trp-Ile EWI motif-containing protein 2; EWI-2; Keratinocytes-associated transmembrane protein 4; KCT-4; LIR-D1; Prostaglandin regulatory-like protein; PGRL; CD316
Human Swiss Prot No.:	Q969P0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD316 Polyclonal Antibody detects endogenous levels of CD316 protein.
Storage and Stability:	-20°C/1 year