

Antibodies Product Datasheet

Product Name:	CD206 Antibody
Catalog Number:	C30572
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human MRC1. AA range:341-390
Concentration:	1 mg/ml
Gene Name:	MRC1
Protein Name:	Macrophage mannose receptor 1
Other Name:	MRC1; CLEC13D; CLEC13DL; MRC1L1; Macrophage mannose receptor 1; MMR; C-type lectin domain family 13 member D; C-type lectin domain family 13 member D-like; Macrophage mannose receptor 1-like protein 1; CD206
Human Swiss Prot No.:	P22897
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD206 Polyclonal Antibody detects endogenous levels of CD206 protein.
Storage and Stability:	-20°C/1 year