

Antibodies Product Datasheet

Product Name:	CD201 Antibody
Catalog Number:	C30574
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Endothelial protein C receptor at AA range: 141-190
Concentration:	1 mg/ml
Gene Name:	PROCR
Protein Name:	Endothelial protein C receptor
Other Name:	PROCR; EPCR; Endothelial protein C receptor; Activated protein C receptor; APC receptor; Endothelial cell protein C receptor; CD201 Q9UNN8
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD201 Polyclonal Antibody detects endogenous levels of CD201 protein.
Storage and Stability:	-20°C/1 year