

Antibodies Product Datasheet

Product Name:	CKR-3 Antibody
Catalog Number:	C30578
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human CCR3. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	CCR3
Protein Name:	C-C chemokine receptor type 3
Other Name:	CCR3; CMKBR3; C-C chemokine receptor type 3; C-C CKR-3; CC-CKR-3; CCR-3; CCR3; CKR3; Eosinophil eotaxin receptor; CD193
Human Swiss Prot No.:	P51677
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CKR-3 Polyclonal Antibody detects endogenous levels of CKR-3 protein.
Storage and Stability:	-20°C/1 year