

Antibodies Product Datasheet

Product Name:	Inhibin β -B Antibody
Catalog Number:	C30587
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Inhibin beta B chain at AA range: 351-400
Concentration:	1 mg/ml
Gene Name:	INHBB
Protein Name:	Inhibin beta B chain
Other Name:	INHBB; Inhibin beta B chain; Activin beta-B chain
Human Swiss Prot No.:	P09529
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Inhibin β -B Polyclonal Antibody detects endogenous levels of Inhibin β -B protein.
Storage and Stability:	-20°C/1 year