

Antibodies Product Datasheet

Product Name:	Inhibin β -A Antibody
Catalog Number:	C30588
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human INHBA. AA range:377-426
Concentration:	1 mg/ml
Gene Name:	INHBA
Protein Name:	Inhibin beta A chain
Other Name:	INHBA; Inhibin beta A chain; Activin beta-A chain; Erythroid differentiation protein; EDF
Human Swiss Prot No.:	P08476
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Inhibin β -A Polyclonal Antibody detects endogenous levels of Inhibin β -A protein.
Storage and Stability:	-20°C/1 year