

Antibodies Product Datasheet

Product Name:	MDA-7 Antibody
Catalog Number:	C30602
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human IL24. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	IL24
Protein Name:	Interleukin-24
Other Name:	IL24; MDA7; ST16; Interleukin-24; IL-24; Melanoma differentiation-associated gene 7 protein; MDA-7; Suppression of tumorigenicity 16 protein
Human Swiss Prot No.:	Q13007
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MDA-7 Polyclonal Antibody detects endogenous levels of MDA-7 protein.
Storage and Stability:	-20°C/1 year