

Antibodies Product Datasheet

Product Name:	IL-29 Antibody
Catalog Number:	C30618
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human IL29. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	IL29
Protein Name:	Interleukin-29
Other Name:	IL29; IFNL1; ZCYTO21; Interleukin-29; IL-29; Cytokine Zcyto21; Interferon lambda-1; IFN-lambda-1
Human Swiss Prot No.:	Q8IU54
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-29 Polyclonal Antibody detects endogenous levels of IL-29 protein.
Storage and Stability:	-20°C/1 year